

Medsys

Gynecology


- HYSTEROSCOPIES



HDScope

- Diameter 2.9 / 4.0 mm
- Viewing Directions 0° / 30°
- Working Length 302 mm
- Standard / Wide Angle
- Compatible to System Storz

- LAPAROSCOPES



HDScope

- Diameter 5.0 / 10.0 mm
- Viewing Directions 0° / 30°
- Working Lengths 312 / 344 mm
- Standard / Wide Angle
- Compatible to System Storz



MEDSYS S.A.
Rue du Gaillard 3
5030 Gembloux

Tel. +32 81 56 84 84
Fax: +32 81 56 88 84
e-mail info@medsys.be
www.medsys.be